

## Why Acero Crowns?

Acero stainless steel resin veneered ES and 100% all zirconia Z1 crowns both exhibit market leading fracture resistance for unbeatable strength while providing the highest esthetic qualities for the patients. Our esthetic line of pediatric crowns simulates the natural translucency and anatomy of the primary dentition resulting in a beautiful restoration.

## Why Zirconia?

Zirconia has been used in dentistry for years and is prized for its strength, biocompatibility, and esthetic qualities.

### Posterior Design

- Flat internal occlusal surface with retentive boxes for easy preparation and strong retention.
- Thin uniform axial thickness for less tooth reduction.
- 24 total anterior crown sizes and 64 total molar sizes to fit any tooth.
- Sandblasted internal surface greatly increases adhesion to the tooth.
- Beautiful esthetics with Vita B1 shade with the translucency and contours of natural teeth.
- User friendly universal numbering system for simple and reliable identification.

### Anterior Design

- Bilateral design allows for less inventory needed.
- Sandblasted internal surface for high adhesion and retention of the crown.
- User friendly Universal Numbering System.
- Beautiful translucency and contours of a natural tooth.

### Acero Zirconia Crowns Details



- Bilateral anterior crowns for less inventory*
- Regular and narrow sizing for molar crowns*
- Universal Numbering System*
- Sandblasted inner surface*
- Uniform axial thickness*
- Feather edge margin*
- Flat occlusal surface*

**ACERO CROWNS** 

[www.acerocrowns.com](http://www.acerocrowns.com)

## ACERO ESTHETIC CROWNS

**ACERO CROWNS**  
[www.acerocrowns.com](http://www.acerocrowns.com)

### ESTHETIC ANTERIOR CROWNS

#### Acero ES and ES Light Esthetic Anterior Pediatric Crowns



### ACERO CROWNS<sup>ES</sup>



1. Acero ES and ES Light are preformed surgical grade stainless steel crowns veneered with a durable composite resin.
2. Acero esthetic crowns come in two shades: Vita B1 for our ES crown and for a whiter shade the ES Light comes in Vita M2 (bleaching shade guide).
3. Acero ES and ESL crowns have been tested by a major US Dental School to have superior fracture, wear, and stain resistance.
4. Acero ES/ESL crowns come in 6 different sizes for the maxillary UR and UL central incisor, lateral incisor and maxillary canine.
5. Like all Acero crowns, they are made from the highest quality surgical grade 316L stainless steel and numbered using the Universal Numbering System A-T.

### ACERO CROWNS<sup>Z1</sup>

#### Acero Z1 Zirconia Crowns

- High Strength
- Load bearing 1200 MPa
- Wear resistant
- Biocompatible
- Superior esthetics



**Developed by Dentists,  
Inspired by Years  
of Clinical Experience**

Trusted, esthetic and affordable dental crowns endorsed by top dentists for making your kid's perfect smile last forever.

### Acero Crowns - ES

SKU	Description
ES D1	Upper right lateral size 1
ES D2	Upper right lateral size 2
ES D3	Upper right lateral size 3
ES D4	Upper right lateral size 4
ES D5	Upper right lateral size 5
ES D6	Upper right lateral size 6
ES E1	Upper right central size 1
ES E2	Upper right central size 2
ES E3	Upper right central size 3
ES E4	Upper right central size 4
ES E5	Upper right central size 5
ES E6	Upper right central size 6
ES F1	Upper left central size 1
ES F2	Upper left central size 2
ES F3	Upper left central size 3
ES F4	Upper left central size 4
ES F5	Upper left central size 5
ES F6	Upper left central size 6
ES G1	Upper left lateral size 1
ES G2	Upper left lateral size 2
ES G3	Upper left lateral size 3
ES G4	Upper left lateral size 4
ES G5	Upper left lateral size 5
ES G6	Upper left lateral size 6
ES U1	Upper canine size 1
ES U2	Upper canine size 2
ES U3	Upper canine size 3
ES U4	Upper canine size 4
ES U5	Upper canine size 5
ES U6	Upper canine size 6
ES KIT	Tray with 1 crown of each size (30 crowns)
ES Tray	Tray only

### Acero Crowns - ES Light

SKU	Description
ESL D1	Upper right lateral size 1
ESL D2	Upper right lateral size 2
ESL D3	Upper right lateral size 3
ESL D4	Upper right lateral size 4
ESL D5	Upper right lateral size 5
ESL D6	Upper right lateral size 6
ESL E1	Upper right central size 1
ESL E2	Upper right central size 2
ESL E3	Upper right central size 3
ESL E4	Upper right central size 4
ESL E5	Upper right central size 5
ESL E6	Upper right central size 6
ESL F1	Upper left central size 1
ESL F2	Upper left central size 2
ESL F3	Upper left central size 3
ESL F4	Upper left central size 4
ESL F5	Upper left central size 5
ESL F6	Upper left central size 6
ESL G1	Upper left lateral size 1
ESL G2	Upper left lateral size 2
ESL G3	Upper left lateral size 3
ESL G4	Upper left lateral size 4
ESL G5	Upper left lateral size 5
ESL G6	Upper left lateral size 6
ESL U1	Upper canine size 1
ESL U2	Upper canine size 2
ESL U3	Upper canine size 3
ESL U4	Upper canine size 4
ESL U5	Upper canine size 5
ESL U6	Upper canine size 6
ESL KIT	Tray with 1 crown of each size (30 crowns)
ESL Tray	Tray only

### Acero Crowns - Z1 Anteriors

SKU	Reference	Description
Z1 DG1	B-U-1	Upper lateral incisor size 1
Z1 DG2	B-U-2	Upper lateral incisor size 2
Z1 DG3	B-U-3	Upper lateral incisor size 3
Z1 DG4	B-U-4	Upper lateral incisor size 4
Z1 DG5	B-U-5	Upper lateral incisor size 5
Z1 DG6	B-U-6	Upper lateral incisor size 6
Z1 EF1	A-U-1	Upper central incisor size 1
Z1 EF2	A-U-2	Upper central incisor size 2
Z1 EF3	A-U-3	Upper central incisor size 3
Z1 EF4	A-U-4	Upper central incisor size 4
Z1 EF5	A-U-5	Upper central incisor size 5
Z1 EF6	A-U-6	Upper central incisor size 6
Z1 CH1	C-U-1	Upper canine size 1
Z1 CH2	C-U-2	Upper canine size 2
Z1 CH3	C-U-3	Upper canine size 3
Z1 CH4	C-U-4	Upper canine size 4
Z1 CH5	C-U-5	Upper canine size 5
Z1 CH6	C-U-6	Upper canine size 6
Z1 MR1	C-L-1	Lower canine size 1
Z1 MR2	C-L-2	Lower canine size 2
Z1 MR3	C-L-3	Lower canine size 3
Z1 MR4	C-L-4	Lower canine size 4
Z1 MR5	C-L-5	Lower canine size 5
Z1 MR6	C-L-6	Lower canine size 6
Z1 Kit Ant		Tray w/2 crowns of ea. size (48 crowns)
Z1 Tray Ant		Replacement tray with no crowns

### Acero Crowns - Z1 Molars

SKU	Reference	Description
Z1 A2	E-UR-2	Upper right primary second molar size 2
Z1 A3	E-UR-3	Upper right primary second molar size 3
Z1 A3N	E-UR-3N	Upper right primary second molar size 3N
Z1 A4	E-UR-4	Upper right primary second molar size 4
Z1 A4N	E-UR-4N	Upper right primary second molar size 4N
Z1 A5	E-UR-5	Upper right primary second molar size 5
Z1 A5N	E-UR-5N	Upper right primary second molar size 5N
Z1 A6	E-UR-6	Upper right primary second molar size 6
Z1 B2	D-UR-2	Upper right primary first molar size 2
Z1 B3	D-UR-3	Upper right primary first molar size 3
Z1 B3N	D-UR-3N	Upper right primary first molar size 3N
Z1 B4	D-UR-4	Upper right primary first molar size 4
Z1 B4N	D-UR-4N	Upper right primary first molar size 4N
Z1 B5	D-UR-5	Upper right primary first molar size 5
Z1 B5N	D-UR-5N	Upper right primary first molar size 5N
Z1 B6	D-UR-6	Upper right primary first molar size 6
Z1 I2	D-UL-2	Upper left primary first molar size 2
Z1 I3	D-UL-3	Upper left primary first molar size 3
Z1 I3N	D-UL-3N	Upper left primary first molar size 3N
Z1 I4	D-UL-4	Upper left primary first molar size 4
Z1 I4N	D-UL-4N	Upper left primary first molar size 4N
Z1 I5	D-UL-5	Upper left primary first molar size 5
Z1 I5N	D-UL-5N	Upper left primary first molar size 5N
Z1 I6	D-UL-6	Upper left primary first molar size 6
Z1 J2	E-UL-2	Upper left primary second molar size 2
Z1 J3	E-UL-3	Upper left primary second molar size 3
Z1 J3N	E-UL-3N	Upper left primary second molar size 3N
Z1 J4	E-UL-4	Upper left primary second molar size 4
Z1 J4N	E-UL-4N	Upper left primary second molar size 4N
Z1 J5	E-UL-5	Upper left primary second molar size 5
Z1 J5N	E-UL-5N	Upper left primary second molar size 5N
Z1 J6	E-UL-6	Upper left primary second molar size 6
Z1 K2	E-LL-2	Lower left primary second molar size 2
Z1 K3	E-LL-3	Lower left primary second molar size 3
Z1 K3N	E-LL-3N	Lower left primary second molar size 3N
Z1 K4	E-LL-4	Lower left primary second molar size 4
Z1 K4N	E-LL-4N	Lower left primary second molar size 4N
Z1 K5	E-LL-5	Lower left primary second molar size 5
Z1 K5N	E-LL-5N	Lower left primary second molar size 5N
Z1 K6	E-LL-6	Lower left primary second molar size 6
Z1 L2	D-LL-2	Lower left primary first molar size 2
Z1 L3	D-LL-3	Lower left primary first molar size 3
Z1 L3N	D-LL-3N	Lower left primary first molar size 3N
Z1 L4	D-LL-4	Lower left primary first molar size 4
Z1 L4N	D-LL-4N	Lower left primary first molar size 4N
Z1 L5	D-LL-5	Lower left primary first molar size 5
Z1 L5N	D-LL-5N	Lower left primary first molar size 5N
Z1 L6	D-LL-6	Lower left primary first molar size 6

### Z1 Molars (Continued)

SKU	Reference	Description
Z1 S2	D-LR-2	Lower right primary first molar size 2
Z1 S3	D-LR-3	Lower right primary first molar size 3
Z1 S3N	D-LR-3N	Lower right primary first molar size 3N
Z1 S4	D-LR-4	Lower right primary first molar size 4
Z1 S4N	D-LR-4N	Lower right primary first molar size 4N
Z1 S5	D-LR-5	Lower right primary first molar size 5
Z1 S5N	D-LR-5N	Lower right primary first molar size 5N
Z1 S6	D-LR-6	Lower right primary first molar size 6
Z1 T2	E-LR-2	Lower right primary second molar size 2
Z1 T3	E-LR-3	Lower right primary second molar size 3
Z1 T3N	E-LR-3N	Lower right primary second molar size 3N
Z1 T4	E-LR-4	Lower right primary second molar size 4
Z1 T4N	E-LR-4N	Lower right primary second molar size 4N
Z1 T5	E-LR-5	Lower right primary second molar size 5
Z1 T5N	E-LR-5N	Lower right primary second molar size 5N
Z1 T6	E-LR-6	Lower right primary second molar size 6
Z1 Kit Molar		Tray with 1 crown of each size (64 crowns)
Z1 Tray Molar		Replacement tray with no crowns

## ES/ES Light & Z1 Anteriors Crowns

### ES & ES Light

Acero esthetic resin veneered stainless steel crowns come in two shades (ES is a Vita B1 shade and the ES Light is similar to the Vita bleaching shade M2).

### Z1 Zirconia

The Z1 is a preformed zirconia crown. Our molar crowns have a flat occlusal inner surface with retentive boxes making tooth preparation simple and quick while providing great retention for the crowns.